



BRED Life Science

## SemenAssay<sup>®</sup> NAG - Kit for Determination of the Neutral alpha-glucosidase Level in Seminal Plasma (Modified Cooper's Method)

The kit is mainly used for evaluating the secretory function of male epididymis, and for topical diagnosis of obstructive azoospermia.

### Features:

- An improved method based on WHO technology
- More convenient and rapider in test with stable results comparing to traditional Cooper method.
- Operation in the microplate. Fast and easy-to-use.
- Available for automatic and semiautomatic tests, good for standardization of test.
- It is ready for use and easy to handle.
- 40 tests per kit and meet the need of batch testing.
- CE approved

